



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : n'Propyl Red
Catalog# : S10157
CAS : 2641-01-2
CI number : -
Formula : $C_{19}H_{23}N_3O_2$
MW : 325.41

TEST	SPECIFICATION
Description	Reddish crystalline powder
Identification (IR)	Conforms to standard
Assay (HPLC)	$\geq 98 \%$
Melting Point	173.5 °C to 175 °C
Visual Transition Interval	From pH 4.6 (red) To pH 6.6 (yellow)
Clarity of Solution: (0.1 g/100 ml – Ethanol)	Clear, red solution

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

