



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Reactive Blue 4
Catalog# : S10165
CAS : 13324-20-4
CI number : 61205
Formula : $C_{23}H_{14}C_{12}N_{10}O_8S_2$
MW : 637.44

TEST	SPECIFICATION
Appearance	Bright Bluish Fine Powder
Dye Content	Between 18 % to 25 %
Lambda Max	595nm (in water)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

