



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Saffron
Catalog# : S10173
CAS : 328.39
CI number : 75100
Formula : $C_{20}H_{24}O_4$
MW : 328.39

TEST	SPECIFICATION
Appearance	Red and Yellow Fibers
Lambda max	427 (452) nm in Ethanol
Solubility	To pass test

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

