



TriStains

**DAWN SCIENTIFIC**

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

## Product Specification

Product : Stains-All  
Catalog# : S10177  
CAS : 7423-31-6  
CI number : -  
Formula :  $C_{30}H_{27}BrN_2S_2$   
MW : 559.60

TEST	SPECIFICATION
Description	Dark-grey powder, Light sensitive
Identification	To pass test
Assay	$\geq 95\%$
Lambda Max	575nm in Ethanol
Solubility	Dark blue to dark purple Solution at 40mg plus 4.0 ml Of Chloroform

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

