



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Sudan Orange G
Catalog# : S10182
CAS : 2051-85-6
CI number : 11920
Formula : $C_{12}H_{10}N_2O_2$
MW : 214.22

TEST	SPECIFICATION
Description	Orange to Red-Orange Powder
Solubility	Very slightly soluble in water
Melting Point	143 - 146 °C
Lambda max (C= 0.008 g/L Methanol)	380-420 in Methanol

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

