



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Sudan Red 7B
Catalog# : S10183
CAS : 6368-72-5
CI number : 26050
Formula : $C_{24}H_{21}N_5$
MW : 379.47

TEST	SPECIFICATION
Appearance	Dark-Red Powder
Identification	To conform to standard
Lambda Max	533 (364) nm in Toluene
Melting Point	130 deg. C
Solubility	Insoluble in water. Slightly soluble in ethanol and soluble in acetone and benzene

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

