



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Victoria Blue B
Catalog# : S10194
CAS : 2580-56-5
CI number : 44045
Formula : $C_{33}H_{32}ClN_3$
MW : 506.10

TEST	SPECIFICATION
Appearance	Dark Blue to Brown powder
Lambda Max	599-605 nm in Methanol
Loss on drying (60 deg. C)	$\leq 5 \%$
Solubility	Slightly soluble in water and somewhat more soluble in the lower aliphatic alcohols

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

