



TriStains

**DAWN SCIENTIFIC**

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

## Product Specification

Product : Victoria Blue R  
Catalog# : S10195  
CAS : 2185-86-6  
CI number : 44040  
Formula :  $C_{29}H_{32}N_3Cl$   
MW : 458.04

TEST	SPECIFICATION
Appearance:	Dark-Blue Powder
Lambda max	615 nm (In Methanol)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

