



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Water Blue 1
Catalog# : S10196
CAS : 28631-66-5
CI number : 42755
Formula : $C_{32}H_{25}N_3Na_2O_9S_3$
MW : 737.73

TEST	SPECIFICATION
Description	Brownish-purple crystalline powder With metallic lustre
Solubility (in water)	Soluble
Lambda max	594 – 610 nm in water

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

