



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Xylidine Ponceau 2R
Catalog# : S10200
CAS : 3761-53-3
CI number : 16150
Formula : $C_{18}H_{14}N_2Na_2O_7S_2$
MW : 480.4

TEST	SPECIFICATION
Appearance	Red Powder
Lambda max	503nm in water
Dye content	~70%
Solubility: 1mg/ml, H ₂ O	Clear red solution

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

